
Software features

Presets	Store and load settings and scale for different objective lens magnifications and stains
Annotations	Various measurements (Arrow, Text, Line, Polyline, Circle, ...)
Auto Exposure	Semi-automatic, automatic
Auto Whitebalance	Semi-automatic, automatic
Snapshot feature	Shortcut to automatically save still images with or without scale bar. Can also be triggered via USB foot switch.
Shading compensation	Active compensation of shading / vignetting
Flicker compensation	Active compensation for older microscopes with flickering light sources

